

Limited of Detection (GT)

Organophosphate

Pesticides	mg./kg.	Pesticides	mg./kg.
acephate	> 30	methamidophos	1
azinphos-ethyl	0.2	methidathion	-ve (test to 150)
azinphos-methyl	> 100	methyl-parathion	> 50
bromophos-rthyl	0.5	monocrotophos	0.05
chlorfenvinphos	0.001	naled (dibrom)	50
chlorpyrifos	0.5	omethoate	20
chlorpyrifos-methyl	> 10	pirimiphos-methyl	10
coumaphos	2	profenophos	0.005
dichlorvos	0.005	prothiophos (tokuthion)	0.1
dimethoate	-ve (test to 150)	phoxim	0.001
diazinon	0.2	phosalone	> 50
dicrotophos	0.05	quinalphos	10
EPN	20 (in methanol)	ronnel	0.01
ethion	0.01	salithion	-ve (test to 20)
fenitrothion	0.5	sulprofos	> 100
isofenphos	> 100	temephos	from 0.01 to 10 samecolor (+1ve)
malathion	8		5
mevinphos	0.005	triazophos	0.05
		trichlorfon	

Carbamate

Pesticides	mg./kg.	Pesticides	mg./kg.
aldicarb	0.1	Carbofuran	0.05
aldicarb-sulphoxide	0.01	Carbosulfan	variable, insoluble in EtOH, have enzyme inhibited.
aldicarb-sulphone	0.2		”
aminocarb	0.1	carbendazim	”
BPMC	5	isoprocarb	0.5
bendiocarb	0.1	methiocarb	> 5
benomyl	variable, insoluble in EtOH, have enzyme inhibited.	methomyl	0.5
		oxamyl	0.05
bufencarb	0.01	pirimicarb	0.1
carbaryl	0.05		



ASIANMEDIC CO., LTD.

Tel: 6689-185-8999, 6690-898-5188

sales@asianmedic.com

www.asianmedic.com