

METHANYL YELLOW TEST KIT




Methyl yellow is a yellow dye for textile, paper, soap and other industries. But in practice, in some developing country these dyes are often misused to color different types of food, among others, noodles, tofu, and yellow snacks. Methyl yellow is irritant so if swallowed can cause irritation of the gastrointestinal tract. The chronic effects of long-term use of methyl yellow on health are liver function, bladder and cancer disorders. Methyl Yellow is a food additives that is prohibited in food use according to Minister of Health Regulation.

PRODUCT INFORMATION

Methyl yellow test kit is a qualitative test kit used to test for methyl yellow dyes prohibited from use as food coloring with food product samples. The principle is

Protect your family from the dangerous of Methanyl Yellow !



**EASY
FAST
ACCURATE
RELIABLE**



PRODUCT INFORMATION

Type	Test Kit
Detection Limit	10 ppm Methyl Yellow
Amount of Test	100 / 200 test
Shelf Life	2 years
Positif Result	The solution color changes to red
Sample	Noodles, tofu and yellow food / food / snacks.
Contents :	<ul style="list-style-type: none">• 1 Reagent MY-1• 1 Reagent MY-2• 3 pcs of pipette

Ordering Information

Cat. No	PAL-016GF (100 test) PAL-017GF (200 test)
---------	--



Asianmedic Co.,Ltd.

Address: 60/3 Soi.Inthamara40,Ratchadaphisak,
Dindeang, Bangkok 10400, Thailand

Hotline: 090.898.5188, 089.185.8999

Tel: 02.691.8348,Fax:02.691.8388

Email: foodtest@asianmedic.com,sales@asianmedic.com